



AIR mobil

Mobile Compressed-AIR
Supply Unit

HEYER 
MEDICAL AG

AIR mobil

Mobile Compressed-AIR Supply Unit

The transportable compressed-AIR supply AIR mobil is a compressor unit for efficient and economic use of the inhalation unit FoggyJet®. It is independent from a stationary compressed-AIR system.

The connection of both devices is made by a pneumatic supply line and an electrical control line. With starting FoggyJet® the AIR mobil is put into operation. After the end of inhalation the AIR mobil stops automatically.

The airflow produced by the compressor is sucked in by a dust/bacteria filter at the front of the unit which is exchangeable. An integrated muffler reduces the noise level so that the unit works quiet.



Part-No.: 194-0000

Technical Data

Electrical supply

Mains voltage	230 V AC
Control voltage	24 V AC
Mains frequency	50 Hz
Capacity	140 VA
International protection	IP 33
Class of protection	I

Performance data compressor

Initial pressure	approx. 1,5 bar
Free flow rate	at least. 19,5 l/min

Ambient conditions

Temperature range	+10° C to +35° C
-------------------	------------------

Dimensions and weight

Height x Breadth x Length	170 x 200 x 300 mm
Weight	6,30 kg

CE 0044

HEYER 
MEDICAL AG

CARL-HEYER-STRASSE 1-3
D-56130 BAD EMS
TEL:+49 2603/791-3-FAX:+49 2603/70424
EMAIL: INFO@HEYERMEDICAL.DE
INTERNET: WWW.HEYERMEDICAL.DE