

**New**  
**New**



# CUMULUS

Modern  
Ultrasonic Nebulization

**HEYER**   
MEDICAL AG

**Wide Range of Application**

HEYER CUMULUS enables a therapy of different infections in almost all areas of the respiratory system. The spectrum reaches from the treatment of the upper airways up to the treatment of the lower respiratory system. Our wide range of accessories guarantees an optimal application at home and in clinics.

**High hygienically comfort**

The cleaning of HEYER Cumulus is very easy. Within seconds the application parts and ultrasonic cup can be dismantled. All parts which are subject to contamination are autoclavable.



0044

**Order-Data**

Part-No.	Description
199-0000	CUMULUS, basic unit, cpl. Heating hose, PPS
199-0200	CUMULUS, basic unit, cpl. Without applications, for sterile-water-system
199-0300	CUMULUS, basic unit, cpl. Infusion-bottel-system INFUDAP
199-0400	CUMULUS, basic unit, cpl. Medicament-system
199-0500	CUMULUS, basic unit, cpl. PPS
100-0100	Ultrasonic cup, autoclavable
500-1097	Permanent post-running System (PPS)
531-2990	Connecting hose, autoclavable
531-2995	Air hose, autoclavable
100-0250	Medicament-system, autoclavable
014-2930	Medicament cup, su = 10 pcs.
198-2001	Heating hose, autoclavable
650-0041	Mouth piece
650-0042	Nose piece
650-0043	Mask, plastics
045-3950	HEYER Cleaning Lotion 1 l for inhalation systems
500-1001	Holding device for patient hose
529-3550	Stand, mobile

Your representative: