

maris^{III}

WELLNESS FOR THE AIRWAYS



WELLNESS FOR THE AIRWAYS

- ❑ Universal inhalation system for the nebulization of brine, thermal and seawater and all kinds of drugs with an indication for inhalation therapy.
- ❑ Separate nebulizer for brine and drug nebulization.
- ❑ Highly integrated, electronically controlled and adjustable heating system for the warm-moist inhalation.
- ❑ Exclusive warming of the brine mist.
- ❑ Type of therapy selectable: continuously drug nebulization, dosed drug nebulization and brine nebulization. The brine and drug nebulization can be used individually or in combination with activation of the aerosol vibration.
- ❑ Therapy time adjustable between 1 to 30 minutes.
- ❑ Built-in aerosol vibration module to provide the aerosol particles with a superimposed oscillation of 100 Hz. This causes a deeper penetration and deposition of mist particles inside the upper airways.
- ❑ Remote control for drug nebulizer to provide an exclusively nebulization of expensive drugs during inspiration.
- ❑ Integrated music module with a connection for headphones.
- ❑ Re-breathing protection to prevent contamination of the nebulizer system by the patient.
- ❑ All patient parts are removable and easy to sterilize, disinfect and autoclave (up to 121°C).
- ❑ Autoclavable stainless-steel bacteria filter for disinfection of the necessary compressed air.
- ❑ HEYER maris III can be mounted with flexible support arm on an inhalation basin, mobile on a mobile stand or used as a tabletop unit.

Basic equipment consisting of:

- Basic unit
- Brine nebulizer
- Drug nebulizer
- Patient part for silicone application parts
- Mouth and nose piece (silicone)
- Stainless steel bacteria filter

Available options:

- Headphones with volume control
- Remote control of the dosed application of drugs
- Patient parts for porcelain application parts
- Various mounting option

