

Technical Specifications									
Type	LuvOn (F1/F2/F3) Full Face Mask			LuvOn (N1/N2/N3) Nasal Mask			LuvOn P1 Nasal Pillows Mask	LuvRes (F1NV/F2NV) Non-Vented Mask	
Model	AF-01	AF-02	AF-03	AN-01	AN-02	AN-03	AP-01	AH-01	AH-02
Therapy Pressure Range	3~40cmH ₂ O			4~30cmH ₂ O			4~20cmH ₂ O	0~40cmH ₂ O	
Flow	(At 50 L/min)	≤0.8cmH ₂ O		≤0.8 cmH ₂ O	≤1.0 cmH ₂ O	≤0.8 cmH ₂ O	≤1.0 cmH ₂ O	≤0.8cmH ₂ O	
Resistance	(At 100 L/min)	≤1.5 cmH ₂ O		≤1.5 cmH ₂ O	≤3.5 cmH ₂ O	≤1.5 cmH ₂ O	≤3.5 cmH ₂ O	≤1.5 cmH ₂ O	
Pressure when Anti-Asphyxia Valve opened	≤2cmH ₂ O			—			—	—	
Pressure when Anti-Asphyxia Valve closed	≤2cmH ₂ O			—			—	—	
Expiratory Resistance with Anti-Asphyxia Valve Open (At 50L/min)	≤2cmH ₂ O			—			—	—	
Inspiratory Resistance with Anti-Asphyxia Valve Open (At 50L/min)	≤2cmH ₂ O			—			—	—	
Noise	≤30dB			≤30dB			≤30dB	≤30dB	



Breathing Mask Series
For Comfortable Breathing



Heyer Care Co.,Ltd.

HQ Add: Building No.6, Zhongguancun
Medical Treatment Park, Daxing District,
Beijing, China

Operation Center Add: No. 10, Chaobai Street,
Yanjiao Development Zone, Sanhe City,
Hebei Province, China

marketing@heyarmed.com



Breathing Mask Series. V1.0 202507





Soft



Lightweight



Stable



Comfortable

- Multi-structure detachable design for easy cleaning and replacement.
- The mask's 6-point adjustable headband ensures a comfortable fit.
- The silicone cushion provides a comfortable seal and minimizes air leaks.
 - The soft, adjustable headband reduces forehead pressure.
 - Stable airflow for a comfortable breathing experience.

LuvOn P1 Nasal Pillows Mask



Relaxed and free with a wide field of vision.

Model: AP-01

- Simplistically designed headgear, comfortable and easy to wear and adjust.
- Ergonomically designed nasal pillows.
- Flexible arm strap support provides an unobstructed view.
- "Quick-release" connector facilitates getting out of bed at night.
- Quiet and lightweight.



Vented Mask Series

Lightweight /Soft

- The forehead strap adjusts to your ideal height for a secure and comfortable fit (AN-01/AN-03).
- The flexible headgear design allows for a wide field of vision.
- Ergonomically designed cushion for optimal comfort.
- "Quick-release" connector (e.g. get out of bed at night).
- Quiet and Comfortable.

LuvOn (N1/N2/N3) Nasal Mask



Patented "butterfly forehead" support technology with adjustable height for comfort and stability.

Comfortable for patients who wear glasses. Suitable for reading enthusiasts and other users who need a wider field of vision.



Simplistically designed for comfort / easy to wear

LuvOn (F1/F2/F3) Full Face Mask



Model: AF-01



Model: AF-02



Model: AF-03

- The unique adjustable forehead support strap can be set to a preferred height, making your experience stable and tight (AF-01/AF-03).
- Flexible arm strap support provides wide vision (AF-02).
- Cushion ergonomically designed.
- "Quick-release" connector (e.g. get out of bed at night).
- Minimalist design.



Non-Vented Mask Series



LuvRes F1NV Non-Vented Mask

LuvRes F2NV Non-Vented Mask



Model: AH-01



Model: AH-02

- Unique elbow design, easily fits for patients in multiple scenarios.
- Cushion ergonomically designed.
- "Quick-release" connector (e.g. get out of bed at night).
- Straight forward design with reliable features.
- The original multi-functional elbow (practical patent) can be used with "exhaust valve/exhaust tee/ double circuit/intubation", which can conveniently meet different clinical needs of doctors.