

LuluNeb Mesh Nebulizer

Designed in Germany

LuluNeb Series, like all other devices from Heyer Medical AG, are developed by highly qualified and experienced engineers using the latest development standards. Production strictly complies with rigorous quality standards.



Corporate Headquarters
heyer Medical AG
Hagenower Str. 73, 19061 Schwerin, Germany
Tel: +49 2603 791 0
Fax: +49 2603 791 200
marketing@heyermid.com
www.heyermid.com

LuluNeb Mesh Nebulizer: V1.1 202606



Breathe Easy, Live Free

LuluNeb Mesh Nebulizer

HOME CARE




LuluNeb
Mesh Nebulizer




www.heyermid.com

Ultra-Light & Hands-Free Convenience

- Ultra-light at just 54g-as light as an egg.
- Wearable design frees your hands.



Medical-Grade Fine Mist for Maximum Efficacy



- Mist particles <5 µm penetrate deep into the lungs.
- Low residue ≤0.2mL ensures every drop counts.
- Nebulization rate ≥0.2 mL/min
- 8mL large cup for uninterrupted treatment.

Whisper-Quiet & Gentle for Kids

- Runs quietly at ≤40dB, quieter than a whisper.



Family-Friendly & All-Scenario

- Includes adult/child masks and a mouthpiece.
- Auto shut-off safely prevents dry without water or solution.
- Type-C charging works anytime, anywhere.

LuluNeb · Breathe Easy, Live Free



Ultra-light & portable (54g)



Hands-free wearable design



Medical-grade fine mist targeting lower airways



Low residual volume ≤0.2mL



Ultra-quiet operation ≤40dB no crying kids



All-in-one solution for the whole family



Easy to clean recharge & use safely



3D skin-friendly mask for comfort

LuluNeb Mesh Nebulizer	
Produc Model	NFS-C
Mesh Sheet Material	Stainless Steel
MMAD	3.8 µm
Percentage of Particles < 5 µm	≥ 60%
Maximum Fill Volume	8 mL
Residual Volume	≤ 0.2 mL
Nebulization Rate	≥ 0.2 mL/min
Aerosol Output Rate	0.072 mL/min
Weight	Main Unit: 36g; Medication Cup: 18g
Dimensions	Main Unit: 5.5 cm x 5.5 cm x 1.8 cm Medication Cup: 3.7 cm x 3.7 cm x 5.1 cm

Important Reminder

For use by individuals aged 3 years and above only.
Always use under the guidance of a healthcare professional.

